

LiteGait[®] 4Home


Product Catalog



LiteGait® 4Home

Designed specifically for individual home use while preserving the exceptional benefits found in clinical LiteGait models. Tailored for purchase or rental by a single user, these Home models are meticulously customized to the user's height, ensuring optimal performance. With a familiar yoke structure, overhead straps, and harness design akin to clinical LiteGait models, the LG4H guarantees superior postural control, enhanced balance, and a secure environment for the practice of various activities including gait, balance exercises, and daily functional tasks.

The LG4H facilitates unweighting, lifting users onto a treadmill or providing essential support for walking and other activities, thus promoting safer and more efficient practice sessions. Additional accessories such as FDY and FlexiStraps are available to enhance versatility, catering to users of all functional levels and enabling dynamic activities to meet individual needs. Whether utilized over ground or on a treadmill like our GaitKeeper Mini, the LG4H ensures an adaptable and effective home rehabilitation experience.



All models are backed by



a one year
replacement
warranty



a responsive
Service
department



an extremely active
Education department
conducting affordable
trainings



exceptional Clinical
Support giving you
direct access to our
clinical educators



LiteGait for Home serves as a versatile tool that adapts to the evolving needs of patients, whether they require intensive practice, ongoing assistance, or maintenance of mobility gains. By empowering patients to practice and progress in the comfort of their own homes, LiteGait promotes long-term independence and well-being.



LiteGait for Home (LG4H) serves as a crucial bridge in enhancing motor learning and rehabilitation outcomes, especially for individuals requiring more practice or assistance.

Home models of LiteGait are ideal for those who would benefit from more practice. Amount of practice is a key factor in motor learning, and most people cannot achieve enough practice within a typical outpatient therapy schedule. For those receiving home-based PT services but requiring a lot of assistance, LG4H will help bridge the gap to allow gait practice in the home environment. Those who have “graduated” from regular PT sessions can benefit from LG4H to maintain gains between episodes of PT care, with intermittent consultation to advance their program

BENEFITS

Benefits for the client:

- ▶ Practice functional tasks in a fall-free environment
- ▶ Increase confidence to engage in activities beyond current functional level
- ▶ Supported environment allows for best quality practice for more efficient, effective movements
- ▶ Enhance symmetry of posture and movement for best practice
- ▶ Vary degree of unweighting individually on each side as needed for decreasing pain, body weight, cardiovascular load, or compensating for weakness
- ▶ Control available planes of movement to increase or decrease the challenge
- ▶ Perform a variety of tasks in different conditions and environments (quadruped, seated, standing, walking, running; bedside, kitchen, outside, treadmill, over ground, with assistive device etc.)




*Be Creative and Challenging
Perform balance and
weight shifting exercises*




shown with Q-pads

LG4Home models combine the affordability and versatility of our LiteGait systems.

Available Options for LG4Home Models



Removable FreeDome—allows the user to rotate 360° or can be locked in any position while safely supported in the LiteGait harness. Practice backward and sideways walking on or off the treadmill or activities such as gait, balance, and ADLs while facing outward.




Tablet BiSym – Objective measurement of the amount and distribution of the support provided by LiteGait. The unique design of LiteGait system allows control of support provided to the patient during therapy to be what the patient needs on each side of the body. One sided weakness or tendencies can be easily corrected and with BiSym the possibly uneven support is measured.





Gaiter Stool – attaches directly to the LiteGait base, either in the back or on either side, to provide positioning and ease for the clinician for facilitating movement over ground or over a treadmill. When working overground, the clinician seated on the GaiterStool can propel the LiteGait by propelling the stool with his or her legs. The Gaiter-Stool is highly adjustable for improved clinician body mechanics.

LiteGait 4 Home Model	LG4H 110	LG4H 220	LG4H 330
Max user weight	110	220	330
Overall height	5'9"	7'0"	7'6"
FreeDome	optional	optional	optional
Max user height overground	5'6"	6'6"	7'0"
Max user height over 6" treadmill	5'0"	6'0"	6'6"
Powered stroke adjustment	18"	24"	28"
Manual max height adjustment	12"	12"	12"
Base width - Standard	27"	27"	31"
Base width - Options	27 30	27 31	31 34
Max Treadmill Width	26 29	26 30	30 33
Door width	31 35	32 36	36.5 39.5
Height of base (floor to top of leg)	5.5"	8"	7"
Length of base leg	82	42.5"	45"
Clearance	13"	18"	20"
Max footprint (width x length)	33" x 32.5"	35" x 42.5"	38.5" x 45"
Harnesses included	1	1	1
BiSym tablet	optional	optional	optional
GaiterStool	optional	optional	optional
AutoStep	optional	optional	optional
Parallel grab bar	included	included	included
Q-Straps	included	included	Included
Power requirements	110V-240V 50/60Hz		





AutoStep— facilitates foot clearance and symmetric steps and reduces the burden of gait facilitation for both pediatric and adult patients. Provides a cue in terminal stance to help forward the leg in swings. Adjustable elastic band attached to the patient's feet is stretched as the stance limb moves back, and shortens to bring the unloaded limb forward in swing.





(LG4H 330)

Includes


- 1 Adult Harness
- Parallel Grab Bar
- Q-Straps



(LG4H 220)

Includes


- 1 Small Harness
- Parallel Grab Bar
- Q-Straps



(LG4H 110)

Includes

- 1 Peds Harness
- Parallel Grab Bar
- Q-Straps




Size chart reflects height without FreeDome yoke

if you fall between sizes, we recommend moving up to the larger model

*If you fall outside of these measurements
please contact a product advisor for more options*

Cost Effective for the Home User

Our traditional LiteGait incorporates essential features needed for a facility to treat large ranges of patient height, weight, and harness girth. These features however, are not needed for a single home user. As a solution, we developed three LiteGait models: LG4H110, LG4H220, and LG4H330. While LG4H units preserve all the design features needed for effective supported therapy, they are customizable for individual use.



*Superior Therapy
& Versatility*

mobility research



480-716-9648



www.LiteGait.com



sales@LiteGait.com